

SpartanTM

**On-site 45 min Legionella
Detection System:
Scientific results**





Mar 22, 2018

Validation overview

- **Objectives**

1. Assess correlation between Spartan qPCR and laboratory culture
2. Compare accuracy of Spartan qPCR vs. laboratory qPCR

- **Categorization of test results**

– <10 GU/mL		Negative
– 10-100 GU/mL		Positive
– 101-1,000 GU/mL		Positive
– >1,000 GU/mL		Positive

Study design

- **51 cooling towers tested for *Legionella pneumophila* for 12 weeks**
 - Spartan qPCR: every 1 week
 - Regular laboratory culture: every 4 weeks
 - Locations: Ottawa, Toronto, Montreal Canada
- **Spartan qPCR**
 - Water samples tested immediately on-site
 - Results in 45 minutes
 - Genomic Units/mL (GU/mL)
- **Laboratory culture**
 - Water samples shipped to laboratory (takes 1-3 days)
 - Results in 10-14 days
 - Colony Forming Units/mL (CFU/mL)






Summary of results

Summary of Spartan qPCR results

- **619 total Spartan qPCR tests**

- 51 towers X 12 weeks X ~1 test/week

- **Spartan qPCR test results (Genomic Units/mL)**

– No result:	44	(7.1%)	
– Undetectable:	444	(72%)	
– <10 GU/mL:	52	(8.4%)	
– 10-100 GU/mL:	65	(11%)	
– 100-1,000 GU/mL :	13	(2.1%)	
– >1,000 GU/mL:	1	(0.16%)	

Objective 1: Spartan qPCR vs. laboratory culture

Spartan qPCR vs. laboratory culture

- **262 water samples tested with both Spartan qPCR and lab culture**
 - Spartan qPCR = Lab culture: 84% concordance
 - Spartan qPCR ≠ Lab culture: 16% discordance
- **Discordant results: why is lab culture under-diagnosing Legionella?**
 - Spartan positive and culture negative: 40/262 (15%)
 - Of the positive results, 40/64 (63%) were missed by culture

	Spartan qPCR (>10 GU/mL)	Spartan qPCR (<10 GU/mL)
Lab culture (>10 CFU/mL)	21	3
Lab culture (<10 CFU/mL)	40	198

Hypothesis: shipping delay causes bacterial growth and degradation

Summary of shipping delay experiments

Relative Change Over Time	Delayed Laboratory Culture	Delayed Spartan qPCR	Delayed Laboratory qPCR
Degradation	67%	72%	77%
Unchanged	22%	16%	9%
Growth	12%	13%	14%

- **Conclusion:** discordant results between Spartan qPCR and laboratory culture are explained by **bacterial growth or degradation** during shipping

Contrived shipping delay

8 HVACs over 3 days
spiked with live
Legionella

Relative Change Over Time	Spiked HVACs, direct Spartan qPCR
Degradation (<2 fold)	66%
No change	23%
Growth (>2 fold)	11%

- Conclusion: Not only are the bacteria dying the DNA is being rapidly degraded in 1-3 days

Objective 2: Spartan qPCR vs. laboratory qPCR

Spartan qPCR vs. laboratory qPCR

- **Spartan qPCR vs. laboratory qPCR**

- Spartan qPCR: water sample tested immediately on site, results in 45 minutes
- Laboratory qPCR: 1-2 day shipping time delay, results in 1 day

- **45 water samples tested with both Spartan qPCR and lab qPCR**

- Spartan qPCR = Lab qPCR: 36% concordance
- Spartan qPCR \neq Lab qPCR: 64% discordance

- **Discordant results: why is lab qPCR under-diagnosing Legionella?**

- Reasons: bacterial degradation due to shipping delay

	Spartan qPCR (> 10 GU/mL)	Spartan qPCR (<10 GU/mL)
Lab qPCR (> 10 GU/mL)	6	3
Lab qPCR (<10 GU/mL)	26	10

Laboratory qPCR vs. laboratory culture

- **Laboratory qPCR vs. laboratory culture**
 - Laboratory qPCR: 1-2 day shipping time delay, results in 1 day
 - Laboratory culture: 1-2 day shipping time delay, results in 10-14 days
- **43 water samples tested with both lab qPCR and lab culture**
 - Lab qPCR = Lab culture: 56% concordance
 - Lab qPCR ≠ Lab culture: 44% discordance
- **Discordant results: lab qPCR is under-diagnosing Legionella**

	Lab qPCR (> 10 GU/mL)	Lab qPCR (< 10 GU/mL)
Lab culture (> 10 CFU/mL)	3	12
Lab culture (< 10 CFU/mL)	7	21

Case study: Cooling tower with Legionella levels >1,000 GU/mL

Case study: culture underdiagnoses Legionella

- **Cooling tower had levels of Legionella >1,000 GU/mL**
 - Lab qPCR and Lab #1 culture: under-diagnosed Legionella
 - Spartan qPCR and Lab #2 culture: diagnosed Legionella


Test type	Week 1	Week 2	Week 3	Week 4	Week 5
Spartan qPCR (GU/mL)	1,300 [†]	980	23	280	240
Direct spike qPCR (GU/mL)	3,100/3,300 [†]	1,100 [‡] /1,700 [*]	-	730 [*]	-
Lab qPCR (GU/mL)	<0.5 [‡] , 8,100	-	-	-	-
Lab #1 culture (CFU/mL)	5 [†]	-	<1	-	-
Lab #2 culture (CFU/mL)	960 [*]	320 [*]	<1 [*]	FedEx Lost	140 [‡]
Spartan Culture (CFU/mL)	11,000 [*]	2,000 [‡]	<4 [‡]	-	-

[†] Time delay of 1 day

[‡] Time delay of 2 days

^{*} Time delay of 3 days

Spartan qPCR: confirmed with direct plating




- **Direct plating**

- HVAC water was directly plated without filter concentration or shipping
- Colonies were counted manually
- Legionella result: 11,000 CFU/mL

- **Confirmatory colony PCR**

- 20 colonies were randomly selected
- PCR performed for 3 targets: *Legionella* species, *L. pneumophila*, *L. pneumophila* serogroup 1
- All 20 colonies were positive for all 3 targets



Why is culture inconsistent: shipping delay + culture method

Different culture methods give different results

- 19 Spartan qPCR positive samples sent for culture
- Samples sent in parallel to 2 different laboratories
- Lab #1 followed the ISO method and the Lab # 2 followed the CDC method
 - Lab #1 was 37X higher on average
- Question: Why the difference between laboratories?

Lab #1 - CDC method (CFU/mL)	Lab #2 - ISO method (CFU/mL)	Lab1/Lab2 Fold difference
520	94	6
500	23	22
320	11	29
120	6	20
80	1	80
60	<1	60
40	<1	40
40	1	40
20	1	20
<1	73	0.01
<1	7	0.1
<1	2	0.5
<1	2	0.5
<1	1	1
<1	<1	1
<1	<1	1
<1	<1	1
<1	<1	1
<1	<1	1
<1	<1	1


Pre-treatments correlate with lower results

- Labs treat samples with acid or heat to remove other bacteria
- ISO methods use pretreatments more often than CDC methods
- Pretreatments frequently associated with lower levels of Legionella
- Conclusion: Culture methods under-report levels of live Legionella

Lab Culture Result (CFU/mL)	Treatment
960	-
520	-
500	-
320	-
120	-
94	Acid
88	Heat
80	-
73	Heat
60	-
40	-
40	-
40	-
23	Acid
20	-
11	Heat
9	Heat
7	Acid
7	Heat
6	Acid
5	-
5	Heat
3	Heat
2	-
2	Acid
2	Heat
1	Acid
1	Acid
1	Acid





Spartan qPCR and direct culture are equivalent

- **Spartan qPCR was directly compared to culture**
 - Samples were tested immediately (no shipping delay)
 - Samples had no pretreatments







- **Conclusion: On-site Spartan qPCR is accurate**

Spartan qPCR and culture correlated well



Lab culture
 Spartan qPCR

Spartan qPCR and culture correlated poorly




Case study of poor correlation

Sampling Date	Spartan qPCR GU/mL	Lab Culture CFU/mL
15-Jan	87	-
12-Jan	36	-
11-Jan	36	-
04-Jan	98	5
03-Jan	250	-
28-Dec	120	-
28-Dec	110	-
21-Dec	63	-
14-Dec	17	<5
07-Dec	76	5
29-Nov	49	<5
22-Nov	40	-
16-Nov	25	-
08-Nov	24	<5


Live Legionella in the water, biofilms and sludge

Source	Spartan qPCR (GU/mL)	Direct Culture (CFU/mL)
Sample point	32	15
Reservoir	31	10
Exit pipe	<LOD	<LOD
Biofilm-Reservoir	Not tested	26 000*
Sludge-Reservoir	1700	3000

* CFU/g as the biofilm was a solid



Biofilm



Sludge

Key findings

- Spartan qPCR results were falsely identified as negative by culture due to bacterial degradation during shipping and culture methods
- Spartan qPCR is equivalent to culture when shipping delay and culture pretreatments are removed
- Case study: one tower had Legionella >1,000 GU/mL that was missed by regularly-scheduled culture and dipslide testing
- On-site culture confirmed high amounts of Legionella in biofilms missed by laboratory culture
- Spartan qPCR is a robust means of monitoring Legionella levels

Main recommendations

- Use on-site qPCR to monitor Legionella levels
- Set 4 action levels for qPCR testing:

Level (GU/mL)	Action
<10	Continue Operation
10-100	Review and Adjust O&M and Water Treatment Program
101-1,000	Clean and Disinfect System
>1,000	Immediately Implement Measures that will Eliminate Water Dispersion by Aerosols